BROOKFIELD

TEXTURE APPLICATION NOTE PEARS

APPLICATION

Simple penetration test of pears.

TEST OBJECTIVE

To quantify the difference between ripe and unripe pears through implementation of simple penetration principles.

BACKGROUND

Pressure tests have traditionally been applied in the quantification of textural properties of pears. Correlation of instrumental mechanical measures with perceived sensory characteristics provides an objective method of quantifying subjective parameters. Quantification of textural properties is beneficial in:

- 1. Determining optimum point of harvest
- 2. Controlled ripening and storage
- 3. Optimization of process operations (e.g., freezing, blanching, cooking, etc.)
- 4. Indicator of potential mechanical injury (bruising)
- 5. Control of raw materials through application of objective quality measures
- 6. Conformance to customer requirements and supply of correctly ripened products
- 7. Quantification of product shelf-life in relation to ripening.

METHOD

Samples were taken from modified atmosphere packaging on day of purchase at ambient conditions. A small section of each fruit was removed (as in illustration) to give stable test surface through which to penetrate.



TABLE 1 LFRA Settings

MODE:	Measure force in	
	compression	
PLOT:	Peak	
SPEED:	0.5mm/s	
DISTANCE:	10	
OPTION:	Normal	
TRIGGER:	Auto 4g Trigger	
	0 00	

PROBE REF: 2mm Ø stainless steel probe, TA 39



SETTING THE STANDARDS in Texture Testing

BROOKFIELD

TEXTURE APPLICATION NOTE PEARS



PARAMETERS HARDNESS

READING

SKIN MODULUS

WORK TO RUPTURE SKIN MEAN OF INTERNAL PEAKS

TOTAL POSITIVE AREA

ADHESIVENESS

ADHESIVE FORCE

Peak positive value within first compression cycle. Direct indication of skin strength. Skin yields and ruptures prior to inner flesh penetration. Curve slope or gradient during penetration. Indicative of flesh rigidity and cell turgidity.

Area of graph.

Mean of positive peaks after skin rupture; indicative of internal cell structure.

Work to complete penetration cycle, including that to initially rupture skin (e.g., sample consistency).

Total negative area, indicative of work required to pull the probe from the sample.

Peak negative force caused by "pinching" and suction of fleshy sample.

RESULTS

	RIPE	UNRIPE
HARDNESS (SKIN) (g)	306	1045
SKIN MODULUS (g/s)	257.1	641.3
WORK TO RUPTURE SKIN (gs)	186.2	866.2
MEAN OF INTERNAL PEAKS (g)	220.8	573.2
TOTAL POSITIVE AREA (gs)	3748.6	11072.2
ADHESIVENESS (gs)	-121.8	-386.1
ADHESIVE FORCE (g)	-87	-327

DISCUSSION

CONCLUSION

The unripe pear is much harder than the ripe pear; it exhibits greater resistance to penetration by the 2mm \emptyset probe. The inner flesh of the ripe pear is much softer than the unripe pear, as shown by the mean of selected peaks, after rupture or yield of the skin.\

This simple test clearly shows differences in pear ripeness.



Brookfield Engineering Laboratories, Inc. 11 Commerce Boulevard, Middleboro, MA 02346-1031 USA Tel: 800-628-8139 or 508-946-6200 • Fax: 508-946-6262 Email: sales@brookfieldengineering.com • Website: www.TextureAnalysis.com