


# Scan® 1200 *High Resolution*

Automatic color colony counter

*Counts up to 7 colors  
on the same plate !*




***Instant, complete & traceable results on all media  
(Petri dishes, Petrifilm™, RIDA™Count / Sanita-kun™, CompactDry™...)***




# Scan<sup>®</sup> 1200: User friendly, simple & high performance

- 1** **Visualisation window**  
Detected colonies are highlighted and marked with a cross
- 2** **Parameter window**  
Only one cursor to optimize counting.  
Customizable settings adaptable to any kind of agar
- 3** **Results window**  
All previously recorded data can be checked again



## Scan<sup>®</sup> 1200 counts in a second



## A full range of accessories



- One cursor to optimize counting
- 6 combinations of lighting and background color; Lighting (upper / lower) / Background (black / white)
- Detects and counts up to 7 colors on the same plate
- Polygonal exclusion areas
- Minimum size of detected colonies: 0,05 mm

- Firewire CCD camera (1280 x 960 pixels)
- Lighting by lifetime lasting LEDs
- English/French/German/Japanese/Chinese/Russian software
- All stainless steel
- Dim : 27 x 27 x 38 cm - 8,7 Kg
- Power : 220 / 110 V - 50 / 60 Hz

Delivered with: Scan<sup>®</sup> software , Firewire PCI card - Option: Firewire PCMCIA card  
Minimum requirements: Windows<sup>™</sup> XP or Vista<sup>™</sup> with Firewire connector (or free PCI slot)

ISO 9001 : 2000 / RoHS / WEEE

**Please fax this form to interscience or to your local distributor**

Scan<sup>®</sup> 1200 colony counter

Quotation


Demonstration

Literature

*I would like more information on:*

- BagMixer<sup>®</sup> 400 ml blender
- MiniMix<sup>®</sup> 100 ml blender
- JumboMix<sup>®</sup> 3500 ml blender
- BagFilter<sup>®</sup> filter bags
- BagTools<sup>®</sup> accessories
- Gravimat<sup>®</sup> dilutor
- Spiral<sup>®</sup> DS+ plater

Name.....  
Company.....  
Address.....  
Zip/City..... Country.....  
Phone..... Fax.....  
email.....  
web.....



**Vitalis Kft**

H-1122 Budapest  
Városmajor u. 22/a

Telefon: +36 (1) 355 4106  
Telefax: +36 (1) 355 4388  
URL: <http://www.vitaliskft.com>

**interscience**

Worldwide: 30, Ch Bois Arpents, 78860 St Nom FRANCE  
Tel. +33(0)1 34 62 62 61 / Fax. +33(0)1 34 62 43 03  
[info@interscience.fr](mailto:info@interscience.fr) / [www.interscience.fr](http://www.interscience.fr)

USA/Canada: 199 Weymouth St. ROCKLAND, MA 02370 USA  
Tel. +1 781 792 2133 / Fax. +1 781 792 2134  
[info@intersciencelab.com](mailto:info@intersciencelab.com) / [www.intersciencelab.com](http://www.intersciencelab.com)